

*Energy scaling of hollow-fiber  
frequency conversion  
of ultrafast pulses*

*Jeff Wojtkiewicz*  
*Charles G. Durfee*

*Support: CSM, NSF*

*OSA Annual Meeting 2002*



Physics Department

Colorado School of Mines

# Hollow fiber frequency mixing



Physics Department


Colorado School of Mines

## Gas as nonlinear medium:


- low dispersion
- supports high intensity
- third-order mixing

## Use capillary to guide pump beams:

- increase interaction length
- Phase-matching of conversion:  
tune pressure to balance gas/waveguide phase



$$\omega_{\text{signal}} = 2 \times \omega_{\text{pump}} - \omega_{\text{idler}}$$




Durfee *et al*, Optics Letters **22**, 1565 (1997)

# Earlier mixing results:



Physics Department

Colorado School of Mines



- pressure optimum produces lowest-order  $HE_{11}$  UV mode
- XPM can *broaden* bandwidth
- compress to 8fs
- high conversion efficiency (40% from pump)
- output energy in the 1-5 $\mu$ J range

**Energy limitations: ionization, cross phase modulation**



# OP-CPA: optical parametric chirped pulse amplification



Physics Department

Colorado School of Mines

Stretch pulse: same intensity, more energy




# OP-CPA: conversion simulations



Physics Department

Colorado School of Mines

## Saturated conversion without gain-narrowing is possible



Shown here:  
4ps pump, 1.8mJ  
50% conversion to UV

- low gain for seed  
→ low degree of gain narrowing
- wide-band or narrow-band pump
- third-order process  
→ generate new wavelengths
- high output beam quality: single mode  
• but requires high input beam quality for efficient waveguide coupling
- output spectrum is not narrowed

# Comparison of amplifier technologies



Physics Department

Colorado School of Mines


	<u>Laser crystal amplifiers</u>	<u>OP-CPA (crystals)</u>	<u>OP-CPA (hollow fiber)</u>
•energy storage	✓	x	x
•ASE/prepulses	✓	x	x
•gain-narrowing	✓	✓	x
•thermal loading	✓	little	x
•dispersion	✓	little	x
•B-integral	some	little	little
•amp $\lambda$	limited choice	$\lambda_a > 400$	$\lambda_a > 100$
<u>pump req'ts</u>			
•pump $\lambda$	limited choice	$\lambda_p > 355$	$\lambda_p > 200$
•duration range	wide	<1ns	<100ps
•shape	any	flat-top	gauss, flat
•beam quality	gauss-flat	flat	gauss
•alignment	easy	overlap, angle	mode-matching
•bandwidth/phase	any	narrow	narrow or chirped

# Hollow fiber OP-CPA: experiment



Physics Department

Colorado School of Mines



- dispersion compensation: short pulse IR into BBO and into capillary
- additional dispersion to blue
- control duration of IR w/add'l material

# Energy extraction




Physics Department

Colorado School of Mines

Without saturation, UV output should *decrease* with greater IR duration:

- signal is linear in IR



With saturation, UV yield increases with IR duration:

- longer IR seed can sweep out more of pump pulse energy
- cross-phase modulation is greatly reduced




# Energy conversion: experimental results



Physics Department

Colorado School of Mines



Input: blue 110 fs, 48μJ      IR: 40μJ

## Output increases with IR chirp:

- better overlap/extraction
- less XPM
- measure 20μJ depletion of blue (may not be seeing all of UV)
- limited by focusability of blue at high conversion

# Dispersion-free SD-FROG





Physics Department

Colorado School of Mines

## Three-mirror self-diffraction FROG:

- simple setup/alignment
- only reflections: wide bandwidth/dispersion-free
- less than 1  $\mu\text{J}$  input required



400nm pulse chirped w/6mm BK7

Durfee *et al*, Optics Letters **24**, 697 (1999)

# High-power doubling




Physics Department

Colorado School of Mines

## Doubling at high conversion efficiency is easy...

- even with “thick” crystals, conversion gives short pulse
- conversion saturates in first layer of crystal



## ... but preserving beam quality is not.


- At high conversion (~50%) focusing quality deteriorates

# Further scaling: opposite chirp mixing



Physics Department

Colorado School of Mines



**Opposite-chirp mixing for blue → ~50% more output (14.5μJ)**

# OP-CPA: compression simulations




Physics Department


Colorado School of Mines

UV output has negative chirp  
(even orders of IR phase change sign)

Compressed UV output as short  
as original IR pulse  
(chirped pump gives shorter output)



Prism pair: +ve dispersion mode  
(correct sign of 2nd and 3rd order phase)



# Future prospects: energy/wavelength scaling




Physics Department

Colorado School of Mines


## Scaling to higher energy:

- higher-energy amplifier
- improve focusability of blue

• current:



• scale to TW level:



## Scaling to short wavelength for high-power VUV:

- pump pulse can be narrow-band
- mix 800nm with: YAG harmonics:
- 355nm → 230nm, 266nm → 160nm, 210nm → 120nm

# Applications of high-power UV pulses



Physics Department

Colorado School of Mines

- pump for recombination XRL
- tunneling ionization in high-frequency limit
- micromachining - materials processing
- hard x-ray generation
- efficient low-order HHG
- photoelectron spectroscopy: fast dynamics of small molecules

# Mixing simulations



Physics Department


Colorado School of Mines

## Propagation code calculates saturated conversion

### Input fields:

- energy, pulse duration, chirp, relative delay

Pressure loop



### Time domain:

- spm, xpm
- nonlinear mixing

### Frequency domain:

- dispersion, losses

### Output processing:

- energy calculation
- post-compression

### Major assumptions:

- *discrete-mode propagation*
- *five harmonic fields*
- *no bending losses*
- *no ionization*



# Characterization of input pulses





Physics Department

Colorado School of Mines

## Blue pulse:

- initially 35fs
- chirped to 110fs with 8.5mm BK7, 3mm fused silica



## IR pulse:

- initially 45 fs
- chirped to 51, 57, 63fs with additional BK7

# Compression




Physics Department

Colorado School of Mines

## UV pulse compression:

- positive IR chirp gives negative UV chirp
- simple compression with material not sufficient: 3rd order
- even orders of phase change sign, odd remain same



- prisms in positive dispersion allow 2nd and 3rd order compensation
- characterize with SD FROG

Group delay for FS prism pair (10cm sep)

